en: INSTRUCTIONS FOR USE	- 2 -
de: GEBRAUCHSANWEISUNG	- 2 -
nl: GEBRUIKSAANWIJZING	- 3 -
fr: MODE D'EMPLOI	- 3 -
es: MODO DE USO	- 4 -
da: BRUGSANVISNING	- 4 -
it: ISTRUZIONI PER L'USO	- 5 -
sv: BRUKSANVISNING	- 5 -
no: BRUKSANVISNING	- 6 -
fi: KÄYTTÖOHJEET	- 6 -
pl: INSTRUKCJA UŻYCIA	- 7 -

[1]

[2]

[3]

[4]

en: INSTRUCTIONS FOR USE

STERILE – do not use if packaging is opened or damaged! For veterinary use only!

System for processing Autologous Conditioned Serum (ACS). Do not use the Orthokine®vet (60 mL) syringe for injections!

1. Blood withdrawal (routine aseptic technique):

Prepare sampling site aseptically. Puncture vein with included butterfly-cannula. Attach Orthokine®vet (60 mL) syringe to butterfly-cannula. Gently draw back plunger of the syringe to the end (click into place) [1]. Wait until blood-flow stops. Detach syringe from butterfly-cannula and close with enclosed green luer lok cap. Gently invert syringe several times. Label Orthokine®vet (60 mL) syringe with patient's name and date and time of blood withdrawal.

2. Incubation:

Incubate blood filled Orthokine®vet (60 mL) syringe at 37° C (98.6° F). Max. time between blood withdrawal and centrifugation: 24 hours!

3. Centrifugation:

Gently invert incubated Orthokine®vet (60 mL) 3 to 5 times before centrifugation. Break off plunger. Place incubated syringe in the centrifuge (load symmetrically with counterbalance). Centrifuge for 10 minutes at 2100 g RCF (RPM according to centrifuge model).

4. Aspiration of ACS (routine aseptic technique):

Gently take Orthokine®vet (60 mL) syringe out of the centrifuge and place in a rack on disinfected working surface. Break off diagonally placed, elongated clips on guide plate [2]. Take off guide plate. Place triangular extraction port over lowest point of serum. Attach 20-30 mL syringe to 100 mm needle. Penetrate the black stopper through triangular extraction port stabilising the white plunger with thumb and index finger [3]. Slowly aspirate ACS. Do not aspirate red blood cells. Gently transfer collected ACS through sterile filter (0.22 µm) and luer luer adapter into aliquot syringes [4]. Close filled aliquot syringes with sterile caps. Label aliquot syringes and freeze (-18° C/-0.4° F or colder) for storage. Shelf life: 7 months post processing. Filter the obtained ACS through sterile filters (0.22 µm) at usage!

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